

Anti-Flotillin 1 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P013140


Product Information

Application	WB, IHC-P/IF (Tissue-P), IP, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Flotillin 1
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Flotillin 1 antibody is for research use only and not for use in diagnostic or therapeutic procedures.


Protein Information

Synonyms	Flotillin-1.
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	O75955
Gene ID	10211
Background	May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Validation Images



Western blot analysis of Flotillin 1 in (1) HeLa lysates; (2) K562 lysates using Flotillin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using Flotillin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.