

Anti-PCSK9 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P011346

Product Information

Application	WB, IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human PCSK9
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PCSK9 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	HCHOLA3, LDLCQ1, NARC1, Neural apoptosis regulated convertase 1, PC9, PCSK9, Proprotein convertase 9.
Calculated MW	Calculated MW: 74 kDa; Observed MW: 62,78 kDa
Uniprot ID	Q8NBP7
Gene ID	255738
Background	May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.



Western blot analysis of PCSK9 in HepG2 lysates using PCSK9 antibody.