

Anti-N WASP Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012960


Product Information

Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA, FC
Reactivity	Human, Mouse
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human WASL
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-N WASP antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	NWASP, WASPB, N-WASP.
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Uniprot ID	O00401
Gene ID	8976
Background	Regulates actin polymerization by stimulating the actin-nucleating activity of the Arp2/3 complex. Binds to HSF1/HSTF1 and forms a complex on heat shock promoter elements (HSE) that negatively regulates HSP90 expression.

Validation Images



Western blot analysis of WASL in HeLa lysates using N WASP antibody.