

## Anti-GABARAPL2 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P013323

### Product Information


Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human GABARAPL2
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-GABARAPL2 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	ATG8, GEF2, ATG8C, GEF-2, GATE16, GATE-16.
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Uniprot ID	P60520
Gene ID	11345
Background	Involved in intra-Golgi traffic. Modulates intra-Golgi transport through coupling between NSF activity and SNAREs activation. It first stimulates the ATPase activity of NSF which in turn stimulates the association with GOSR1.

## Validation Images

---



Western blot analysis of GABARAPL2 in (1) HeLa lysates; (2) mouse spleen lysates; (3) rat brain lysates using GABARAPL2 antibody.