

## Anti-PKC theta Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P010835

### Product Information


Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human PKC theta
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PKC theta antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	PRKCQ, PRKCT, Protein kinase C theta type, nPKC-theta.
Calculated MW	Calculated MW: 82 kDa; Observed MW: 82 kDa
Uniprot ID	Q04759
Gene ID	5588
Background	This is a calcium-independent, phospholipid-dependent, serine- and threonine-specific enzyme. Essential for T-cell receptor (TCR)-mediated T-cell activation, but is dispensable during TCR-dependent thymocyte development.

## Validation Images

---



Western blot analysis of PKC theta in Jurkat lysates using PKC theta antibody.