

Anti-Prostate Secretory Protein Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012961


Product Information

Application	WB, IHC-P/IF (Tissue-P), IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Prostate Secretory Protein/PSP
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Prostate Secretory Protein antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	MSP, PSP, IGBF, MSPB, PN44, PRPS, HPC13, PSP57, PSP94, PSP-94.
Calculated MW	Calculated MW: 13 kDa; Observed MW: 13 kDa
Uniprot ID	P08118
Gene ID	4477
Background	Specific receptors for this protein are found on spermatozoa and in the prostate.

Validation Images



Western blot analysis of Prostate Secretory Protein/PSP in Human prostate cancer lysates using Prostate Secretory Protein antibody.