

## Anti-PAX5 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P012887

### Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Peptide sequence around derived from Human Pax5.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PAX5 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	PAX5, Paired box protein Pax-5, B-cell-specific transcription factor, BSAP.
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	Q02548
Gene ID	5079
Background	PAX5 May play an important role in B-cell differentiation as well as neural development and spermatogenesis. Involved in the regulation of the CD19 gene, a B-lymphoid-specific target gene. Interacts with DAXX. Binds DNA as a monomer. Binds TLE4.



Western blot analysis of PAX5 in MOLT4 and JK lysates using Pax5 antibody.