

## Anti-PKNOX2 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P010865

### Product Information


Application	WB, IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Mouse
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthesized peptide derived from the C-terminal region of human PREP-2.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PKNOX2 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	PKNOX2, PREP2, Homeobox protein PKNOX2, Homeobox protein PREP-2, PBX/knotted homeobox 2.
Calculated MW	Calculated MW: 52 kDa; Observed MW: 55 kDa
Uniprot ID	Q96KN3
Gene ID	63876
Background	Homeodomain proteins are sequence-specific transcription factors that share a highly conserved DNA-binding domain and play fundamental roles in cell proliferation, differentiation, and death.

## Validation Images

---



Western blot analysis of PKNOX2 in HepG2 lysates using PKNOX2 antibody.