

Anti-Cav gamma 5 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012553


Product Information

Application	WB, IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthetic Peptide of L-type Ca ⁺⁺ CP γ 5
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Cav gamma 5 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	acng5 antibody, TARP gamma-5 antibody, Transmembrane AMPAR regulatory protein gamma-5 antibody.
Calculated MW	Calculated MW: 31 kDa; Observed MW: 35 kDa
Uniprot ID	Q9UF02
Gene ID	27091
Background	Regulates the gating properties of AMPA-selective glutamate receptors (AMPA receptors).

Validation Images



Western blot analysis of Cav gamma 5 in Human Brain, mouse Brain, rat Brain lysates using Cav Cav gamma 5 antibody.