

Anti-Phospho-BRAF (Thr401) Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P010588


Product Information

Application	WB, ICC/IF (Cell), IP, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IF 1:50~1:200; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Phospho-B Raf (T401)
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-BRAF (Thr401) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	BRAF, BRAF1, RAFB1, Serine/threonine-protein kinase B-raf, Proto-oncogene B-Raf, p94, v-Raf murine sarcoma viral oncogene homolog B1.
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Uniprot ID	P15056
Gene ID	673
Background	Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron. Frequently mutated in thyroid cancers, skin melanomas and at lower frequency in a wide range of human cancers. An activating mutation, mimicking phosphorylation of the activation loop, is seen in 60% of malignant melanoma samples. Raf mutations are generally exclusive to Ras activating mutations.

Validation Images



Western blot analysis of Phospho-B Raf (T401) in (1) PC-12 lysates; (2) PC-12 lysates treated with TPA using Phospho-BRAF (Thr401) antibody.