

Anti-Wnt16 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P010683


Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthesized peptide derived from the Internal region of human Wnt-16.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Wnt16 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	WNT16, Protein Wnt-16.
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Uniprot ID	Q9UBV4
Gene ID	51384
Background	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters .

Validation Images



Western blot analysis of Wnt16 in rat kidney lysates using Wnt16 antibody.