

## Anti-IL-17B Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012742

### Product Information


Application	IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Rat
Dilution	IHC-P 1:50~1:100; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthetic peptide from human protein.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-IL-17B antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	IL-17B, interleukin 17B, NIRF, IL-20, ZCYTO7.
Uniprot ID	Q9UHF5
Gene ID	27190
Background	Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.

## Validation Images

---



Immunohistochemistry analysis of paraffin-embedded Human smallintestine using IL-17B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.