

## Anti-C1RL Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P013359

### Product Information


Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthesized peptide derived from human C1RL.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-C1RL antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	Complement C1r subcomponent-like protein (C1r-LP, C1r-like protein, EC 3.4.21.-, C1r-like serine protease analog protein, CLSPa).
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Uniprot ID	Q9NZP8
Gene ID	51279
Background	Mediates the proteolytic cleavage of HP/haptoglobin in the endoplasmic reticulum.

## Validation Images

---



Western blot analysis of C1RL in HeLa lysates using C1RL antibody.