

# Anti-TNF alpha Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M011300

## Product Information


Application	WB, IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Rat, Mouse (Cell)
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100
Host	Mouse
Clonality	Monoclonal
Clone No.	54K58L38
Isotype	IgG1
Label	Unconjugated
Immunogen	Purified recombinant protein expressed in E.coli.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-TNF alpha antibody [54K58L38] is for research use only and not for use in diagnostic or therapeutic procedures.

## Protein Information


Synonyms	TNF, TNFA, TNFSF2, Tumor necrosis factor, Cachectin, TNF-alpha, Tumor necrosis factor ligand superfamily member 2, TNF-a.
Calculated MW	Calculated MW: 26 kDa; Observed MW: 17,25 kDa
Uniprot ID	P01375
Gene ID	7124
Background	This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR2. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.

## Validation Images

---



Western blot analysis of TNF alpha in Recombinant Human TNF a Protein using TNF alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Tissue using TNF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.