

Anti-Phospho-Rac1 (Ser71) Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P011183


Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthetic peptide of human RAC1
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-Rac1 (Ser71) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	Ras-related C3 botulinum toxin substrate 1, Cell migration-inducing gene 5 protein, Ras-like protein TC25, p21-Rac1, RAC1, TC25.
Calculated MW	Calculated MW: 21 kDa; Observed MW: 28 kDa
Uniprot ID	P63000
Gene ID	5879
Background	Plasma membrane-associated small GTPase which cycles between active GTP-bound and inactive GDP-bound states.

Validation Images



Western blot analysis of Phospho-Rac1 (Ser71) in 293 lysates using Phospho-Rac1 (Ser71) antibody.