

Anti-Bcl10 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R013095


Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	17M59M70
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human Bcl10
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Bcl10 antibody [17M59M70] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	BCL10, CIPER, CLAP, B-cell lymphoma/leukemia 10, B-cell CLL/lymphoma 10, Bcl-10, CARD-containing molecule enhancing NF-kappa-B, CARD-like apoptotic protein, hCLAP, CED-3/ICH-1 prodomain homologous E10-like regulator, CIPER, Cellular homolog.
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Uniprot ID	O95999
Gene ID	8915
Background	Promotes apoptosis, pro-caspase-9 maturation and activation of NF-kappa-B via NIK and IKK. May be an adapter protein between upstream TNFR1-TRADD-RIP complex and the downstream NIK-IKK-IKAP complex. Is a substrate for MALT1.

Validation Images



Western blot analysis of Bcl10 in K562, C6, 3T3, HeLa lysates using Bcl10 antibody.