

Anti-CDC25 Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R011516

Product Information

Application	WB, ELISA
Reactivity	Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	98L88L09
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human CDC25
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-CDC25 antibody [98L88L09] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	GGRP, GRF1, CDC25, GRF55, CDC25L, H-GRF55, PP13187, ras-GRF1.
Calculated MW	Calculated MW: 145 kDa; Observed MW: 145 kDa
Uniprot ID	Q13972
Gene ID	5923
Background	Promotes the exchange of Ras-bound GDP by GTP.



Western blot analysis of CDC25 in rat Brain lysates using CDC25 antibody.