

Anti-FYB Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011167


Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	44M36K26
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human FYB
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-FYB antibody [44M36K26] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	ADAP, Fyb, FYN binding protein, p120/p130, SLAP130.
Calculated MW	Calculated MW: 85 kDa; Observed MW: 130,120 kDa
Uniprot ID	O15117
Gene ID	2533
Background	Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

Validation Images



Western blot analysis of FYB in Jurkat lysates using FYB antibody.