

Anti-GABARAPL2 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011546


Product Information

Application	WB, ELISA
Reactivity	Human, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	97K29M37
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human GABARAPL2
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-GABARAPL2 antibody [97K29M37] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	ATG8, GEF2, ATG8C, GEF-2, GATE16, GATE-16.
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Uniprot ID	P60520
Gene ID	11345
Background	Involved in intra-Golgi traffic. Modulates intra-Golgi transport through coupling between NSF activity and SNAREs activation. It first stimulates the ATPase activity of NSF which in turn stimulates the association with GOSR1.

Validation Images



Western blot analysis of GABARAPL2 in C6, HeLa lysates using GABARAPL2 antibody.