

Anti-ICAD Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011443


Product Information

Application	WB, ELISA
Reactivity	Human, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	15K94L14
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human ICAD
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-ICAD antibody [15K94L14] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	DFFA, DFF1, DFF45, H13, DNA fragmentation factor subunit alpha, DNA fragmentation factor 45 kDa subunit, DFF-45, Inhibitor of CAD, ICAD.
Calculated MW	Calculated MW: 37 kDa; Observed MW: 45 kDa
Uniprot ID	O00273
Gene ID	1676
Background	Inhibitor of the caspase-activated DNase (DFF40).

Validation Images



Western blot analysis of ICAD in rat Brain, MCF-7, C6 lysates using ICAD antibody.