

Anti-PCSK9 Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R012788

Product Information

Application	WB, IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	52K15L97
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human PCSK9
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PCSK9 antibody [52K15L97] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	HCHOLA3, LDLCQ1, NARC1, Neural apoptosis regulated convertase 1, PC9, PCSK9, Proprotein convertase 9.
Calculated MW	Calculated MW: 74 kDa; Observed MW: 65,80 kDa
Uniprot ID	Q8NBP7
Gene ID	255738
Background	May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.



Western blot analysis of PCSK9 in HeLa lysates using PCSK9 antibody.