

Anti-RhoGDI Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R012648


Product Information

Application	WB, ELISA
Reactivity	Mouse, Rat, Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	49L51L33
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human RhoGDI
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-RhoGDI antibody [49L51L33] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	ARHGDI, GDIA1, Rho GDP-dissociation inhibitor 1, Rho GDI 1, Rho-GDI alpha.
Calculated MW	Calculated MW: 23 kDa; Observed MW: 26 kDa
Uniprot ID	P52565
Gene ID	396
Background	Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them.

Validation Images



Western blot analysis of RhoGDI in K562, rat Brain, C6, 3T3, HeLa lysates using RhoGDI antibody.