

Anti-ETS1 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R013077


Product Information

Application	WB, ELISA
Reactivity	Human, Mouse
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	85L02L76
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human ETS1
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-ETS1 antibody [85L02L76] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	ETS1, EWSR2, Protein C-ets-1, p54.
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Uniprot ID	P14921
Gene ID	2113
Background	Transcription factor.

Validation Images



Western blot analysis of ETS1 in HeLa, A549, HL-60, U2OS lysates using ETS1 antibody.