

## Anti-Bcl2 L14 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R013393


### Product Information

Application	WB, ICC/IF (Cell), IP, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IF 1:50~1:200; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	37K68M46
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of human Bcl G/BCL2L14
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Bcl2 L14 antibody [37K68M46] is for research use only and not for use in diagnostic or therapeutic procedures.


### Protein Information

Synonyms	BCLG.
Calculated MW	Calculated MW: 37 kDa; Observed MW: 43 kDa
Uniprot ID	Q9BZR8
Gene ID	79370
Background	Plays a role in apoptosis.


## Validation Images



Western blot analysis of Bcl G/BCL2L14 in C6, 3T3, PC-12, Raw264.7 lysates using Bcl L14 antibody.



Western blot analysis of Bcl G/BCL2L14 in HeLa, A549, HL-60, U251, U87-MG lysates using Bcl G/BCL2L14 antibody.



Immunocytochemistry analysis of Bcl G/BCL2L14 (green) in HeLa cells using Bcl G/BCL2L14 antibody and DAPI (blue)