

Anti-Bora Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R013142


Product Information

Application	WB, IP, ELISA
Reactivity	Human, Rat
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	80L70K82
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human Bora
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Bora antibody [80L70K82] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	C13orf34.
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Uniprot ID	Q6PGQ7
Gene ID	79866
Background	Required for the activation of AURKA at the onset of mitosis.

Validation Images



Western blot analysis of Bora in HL-60, C6 lysates using Bora antibody.