

Anti-CD229 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011294


Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	96M06M37
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of human CD229
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-CD229 antibody [96M06M37] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	hly9, mLY9, CD229, SLAMF3.
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Uniprot ID	Q9HBG7
Gene ID	4063
Background	Self-ligand receptor of the signaling lymphocytic activation molecule (SLAM) family.

Validation Images



Western blot analysis of CD229 in HeLa, A549, HL-60, U251 lysates using CD229 antibody.