

## Anti-HNF 4 alpha Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R013485

### Product Information


Application	WB, IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	71L08K69
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of human HNF-4-alpha
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-HNF 4 alpha antibody [71L08K69] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	TCF, HNF4, MODY, FRTS4, MODY1, NR2A1, TCF14, HNF4a7, HNF4a8, HNF4a9, NR2A21, HNF4alpha.
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Uniprot ID	P41235
Gene ID	3172
Background	HNF-4-alpha is a transcription factor that belongs to the steroid hormone receptor superfamily. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha.

## Validation Images

---



Western blot analysis of HNF4alpha in A549, MCF-7 lysates using HNF 4 alpha antibody.