

## Anti-LILRB1 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R012223

### Product Information


Application	WB, IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	81K04K29
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of human LILRB1
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-LILRB1 antibody [81K04K29] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	ILT2, LIR1, MIR7, PIRB, CD85J, ILT-2, LIR-1, MIR-7, PIR-B.
Calculated MW	Calculated MW: 71 kDa; Observed MW: 110 kDa
Uniprot ID	Q8NHL6
Gene ID	10859
Background	Receptor for class I MHC antigens.

## Validation Images

---



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.