

## Anti-SAP97 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011700

### Product Information


Application	WB, ELISA
Reactivity	Mouse, Rat, Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	73K07L06
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human SAP97
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-SAP97 antibody [73K07L06] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information


Synonyms	DLG1, DLGH1, hDlg, SAP97.
Calculated MW	Calculated MW: 100 kDa; Observed MW: 100 kDa
Uniprot ID	Q12959
Gene ID	1739
Background	Essential multidomain scaffolding protein required for normal development (By similarity). Recruits channels, receptors and signaling molecules to discrete plasma membrane domains in polarized cells.

## Validation Images

---



Western blot analysis of SAP97 in mouse brain, rat brain lysates using SAP97 antibody.



Western blot analysis of SAP97 in HeLa, A549, HL-60 lysates using SAP97 antibody.