

## Anti-DAP kinase 3 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011575

### Product Information


Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	27K90K37
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human ZIP Kinase
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-DAP kinase 3 antibody [27K90K37] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	DAPK3, ZIPK, Death-associated protein kinase 3, DAP kinase 3, DAP-like kinase, Dlk, MYPT1 kinase, Zipper-interacting protein kinase, ZIP-kinase.
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Uniprot ID	O43293
Gene ID	1613
Background	Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.

## Validation Images

---



Western blot analysis of ZIP Kinase in HeLa, A549 lysates using DAP kinase 3 antibody.