

Anti-PIAS1/2 Rabbit mAb

Purified Rabbit Monoclonal Antibody

Catalog # R011803


Product Information

Application	WB, IP, ELISA, IHC-P/IF (Tissue-P)
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	60K48M21
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human PIAS1
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PIAS1/2 antibody [60K48M21] is for research use only and not for use in diagnostic or therapeutic procedures.


Protein Information

Synonyms	PIAS1, PIAS2, Protein inhibitor of activated STAT protein 1.
Calculated MW	Calculated MW: 72 kDa; Observed MW: 76 kDa
Uniprot ID	O75925, O75928
Gene ID	8554/9063
Background	This proteins are involved in the pathway protein sumoylation, which is part of Protein modification.

Validation Images



Western blot analysis of PIAS1/PIAS2 in mouse testis, rat testis lysates using PIAS1/2 antibody.



Western blot analysis of PIAS1+PIAS2 in HeLa, A549, HL-60, U2OS lysates using PIAS1+PIAS2 antibody.