

## Anti-GAB1 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012410

### Product Information


Application	WB, ELISA
Reactivity	Human, Mouse
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Peptide sequence around derived from Human Gab1.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-GAB1 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information


Synonyms	GAB1, GRB2-associated-binding protein 1, GRB2-associated binder 1, Growth factor receptor bound protein 2-associated protein 1.
Calculated MW	Calculated MW: 77 kDa; Observed MW: 110 kDa
Uniprot ID	Q13480
Gene ID	2549
Background	Probably involved in EGF and insulin receptor signaling.

## Validation Images

---



Western blot analysis of GAB1 in A431 and MCF lysates using Gab1 (Ab627) antibody.



Western blot analysis of GAB1 in A431 and MCF-7 lysates using Gab1 (Ab627) antibody.