

## Anti-Cyclin A1 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P011184

### Product Information

Application	WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA
Reactivity	Monkey, Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	The antiserum was produced against synthesized peptide derived from human Cyclin A.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Cyclin A1 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	CCNA1, Cyclin-A1.
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Uniprot ID	P78396
Gene ID	8900
Background	May be involved in the control of the cell cycle at the G1/S (start) and G2/M (mitosis) transitions.



Western blot analysis of Cyclin A1 in various lysates using Cyclin A1 antibody.