

Anti-SREBP1 Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P010134


Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthesized peptide derived from the Internal region of human SREBP-1.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-SREBP1 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	SREBF1, BHLHD1, SREBP1, Sterol regulatory element-binding protein 1, SREBP-1, Class D basic helix-loop-helix protein 1, bHLHD1, Sterol regulatory element-binding transcription factor 1.
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Uniprot ID	P36956
Gene ID	6720
Background	Sterol regulatory element-binding protein 1: Precursor of the transcription factor form (Processed sterol regulatory element-binding protein 1), which is embedded in the endoplasmic reticulum membrane.

Validation Images



Western blot analysis of SREBP1 in HeLa, HepG2 lysates using SREBP1 antibody.