

## Anti-BRAP Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M013029

### Product Information


Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Mouse
Clonality	Monoclonal
Clone No.	72L37L66
Isotype	IgG1
Label	Unconjugated
Immunogen	Purified recombinant human BRAP protein fragments expressed in E.coli.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-BRAP antibody [72L37L66] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information


Synonyms	IMP, BRAP2, RNF52.
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Uniprot ID	Q7Z569
Gene ID	8315
Background	Negatively regulates MAP kinase activation by limiting the formation of Raf/MEK complexes probably by inactivation of the KSR1 scaffold protein.

## Validation Images

---



Western blot analysis of BRAP in CHO-K1 lysates, transfected with a Human pFLAG-CMV-2-BRAP construct (A) or transfected with a Human pFLAG-CMV-2 construct (B), using BRAP antibody [72L37L66].



Western blot analysis of BRAP in A549 and HUVEC lysates using BRAP antibody [72L37L66].