

Anti-Flotillin 1 Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M010587


Product Information

Application	WB, ELISA
Reactivity	Mouse (Cell), Rat
Dilution	WB 1:500~1:1,000
Host	Mouse
Clonality	Monoclonal
Clone No.	33K01M26
Isotype	IgG1
Label	Unconjugated
Immunogen	Synthetic Peptide of Flotillin-1
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Flotillin 1 antibody [33K01M26] is for research use only and not for use in diagnostic or therapeutic procedures.


Protein Information

Synonyms	Flotillin-1.
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	O75955
Gene ID	10211
Background	May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Validation Images



Western blot analysis of Flotillin 1 in rat Brain, mouse Brain lysates using Flotillin1 antibody [33K01M26].



Western blot analysis of Flotillin 1 in mouse brain, rat brain lysates using Flotillin1 antibody [33K01M26].