

Anti-Bcl2 Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M010582


Product Information

Application	WB, ELISA
Reactivity	Human, Rat, Mouse (Cell)
Dilution	WB 1:500~1:1,000
Host	Mouse
Clonality	Monoclonal
Clone No.	26L16K77
Isotype	IgG1
Label	Unconjugated
Immunogen	Synthetic peptide conjugated to KLH.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Bcl2 antibody [26L16K77] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	BCL2, Apoptosis regulator Bcl-2.
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Uniprot ID	P10415
Gene ID	596
Background	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma.

Validation Images



Western blot analysis of BCL2 in THP-1 lysates using BCL2 mouse Monoclonal Antibody [26L16K77].