

Anti-GRB2 Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M010448


Product Information

Application	WB, ELISA
Reactivity	Mouse (Cell), Rat, Human
Dilution	WB 1:500~1:1,000
Host	Mouse
Clonality	Monoclonal
Clone No.	52L53L01
Isotype	IgG1
Label	Unconjugated
Immunogen	Recombinant protein of human GRB2.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-GRB2 antibody [52L53L01] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	GRB2, ASH, Growth factor receptor-bound protein 2, Adapter protein GRB2, Protein Ash, SH2/SH3 adapter GRB2.
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Uniprot ID	P62993
Gene ID	2885
Background	Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

Validation Images



Western blot analysis of GRB2 in K562, Raji, Jurkat, Ramos, C6, 3T3 and HeLa lysates using GRB2 antibody.