

Anti-Lyn Mouse mAb

Purified Mouse Monoclonal Antibody

Catalog # M013380


Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Mouse
Clonality	Monoclonal
Clone No.	27K94M01
Isotype	IgG1
Label	Unconjugated
Immunogen	Purified recombinant human Lyn protein fragments expressed in E.coli.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Lyn antibody [27K94M01] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	LYN, JTK8, Tyrosine-protein kinase Lyn, Lck/Yes-related novel protein tyrosine kinase, V-yes-1 Yamaguchi sarcoma viral related oncogene homolog, p53Lyn, p56Lyn.
Calculated MW	Calculated MW: 59 kDa; Observed MW: 56,53 kDa
Uniprot ID	P07948
Gene ID	4067
Background	Lyn a tyrosine kinase of the Src family that plays a role in B cell antigen-receptor signaling.

Validation Images



Western blot analysis of Lyn in Ramos, Raji, K562 and HeLa lysates using Lyn antibody.