

Anti-mTOR Mouse mAb

Purified Recombinant Mouse Monoclonal Antibody


Catalog # M900010

Product Information

Application	WB, IHC-P/IF (Tissue-P), IF (Cell)/ICC, ELISA
Reactivity	Human, Mouse (Cell), Rat
Dilution	WB 1:1,000~1:5,000; IHC-P 1:100~1:200; IF 1:100~1:200
Host	Mouse
Clonality	Monoclonal
Clone No.	91K70Q97
Isotype	IgG
Label	Unconjugated
Immunogen	This mTOR antibody is generated from mice immunized with a recombinant human mTOR protein expressed in <i>E. coli</i> .
Format	Affinity purified monoclonal antibody supplied in PBS with 0.01% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-mTOR Mouse mAb [91K70Q97] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	dJ576K7.1 (FK506 binding protein 12 rapamycin associated protein 1), FK506 binding protein 12 rapamycin associated protein 1, FK506 binding protein 12 rapamycin associated protein 2, FK506 binding protein 12 rapamycin complex associated protein 1, FK506-binding protein 12-rapamycin complex-associated protein 1, FKBP rapamycin associated protein, FKBP12 rapamycin complex associated protein, FKBP12-rapamycin complex-associated protein 1, FKBP12-rapamycin complex-associated protein, FLJ44809, FRAP, FRAP1, FRAP2, Mammalian target of rapamycin, Mechanistic target of rapamycin, mTOR, MTOR_HUMAN, OTTHUMP0000001983, RAFT1, Rapamycin and FKBP12 target 1, Rapamycin associated protein FRAP2, Rapamycin target protein 1, Rapamycin target protein, RAPT1, Serine/threonine-protein kinase mTOR.
Calculated MW	Calculated MW: 289 kDa; Observed MW: 289 kDa
Uniprot ID	P42345, P42346, Q9JLN9
Gene ID	2475
Background	FRAP1 belongs to a family of phosphatidylinositol kinase-related kinases. These kinases mediate cellular responses to stresses such as DNA damage and nutrient deprivation. This protein acts as the target for the cell-cycle arrest and immunosuppressive effects of the FKBP12-rapamycin complex. FRAP1 is a part of the TORC2 complex which plays a critical role in AKT1 Ser-473 phosphorylation, and may modulate the phosphorylation of PKCA and regulate actin cytoskeleton organization.



Western Blot - Anti-mTOR Mouse mAb [91K70Q97]

All lanes: M900010 at 1:5,000 dilution

Lane 1: HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates

Lane 2: HepG2 (Human hepatocellular carcinoma epithelial cell) whole cell lysates

Lane 3: HCT116 (Human colorectal carcinoma epithelial cell) whole cell lysates

Lane 4: A431 (Human epidermoid carcinoma cell) whole cell lysates

Lane 5: T24 (Human bladder cancer epithelial cell) whole cell lysates

Lane 6: RAW264.7 (Mouse mononuclear macrophage leukemia epithelial cell) whole cell lysates


Lysates/proteins at 10 µg per lane.

Secondary antibody: Goat Anti-Mouse IgG(H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution

Predicted band size: 289 kDa

Observed band size: 289 kDa

Developed using the ECL technique (Cat. No. SQ201).



Immunofluorescence - Anti-mTOR Mouse mAb [91K70Q97]


Sample: HeLa cells

The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.

Primary antibodies: M900010 at 1:100 dilution

Secondary antibodies: Goat anti-Mouse (488) at 1:1,000 dilution (shown in green)

Nuclei were stained with DAPI (shown in blue).



Immunofluorescence - Anti-mTOR Mouse mAb [91K70Q97]


Sample: HeLa cells

The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.

Primary antibodies: M900010 at 1:100 dilution

Secondary antibodies: Goat anti-Mouse (488) at 1:1,000 dilution (shown in green)

Nuclei were stained with DAPI (shown in blue).



Immunofluorescence - Anti-mTOR Mouse mAb [91K70Q97]


Sample: HeLa cells

The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.

Primary antibodies: M900010 at 1:100 dilution

Secondary antibodies: Goat anti-Mouse (488) at 1:1,000 dilution (shown in green)

Nuclei were stained with DAPI (shown in blue).



Immunofluorescence - Anti-mTOR Mouse mAb [91K70Q97]


Sample: HeLa cells

The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.

Primary antibodies: M900010 at 1:100 dilution

Secondary antibodies: Goat anti-Mouse (488) at 1:1,000 dilution (shown in green)

Nuclei were stained with DAPI (shown in blue).



Immunohistochemistry - Anti-mTOR Mouse mAb [91K70Q97]

Sample: Paraformaldehyde-fixed, paraffin embedded rat lung tissue

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.

Primary antibody: M900010 at 1:200 dilution

Secondary antibody: Goat Anti-Mouse IgG (H+L), HRP conjugated at 1:1,000 dilution

DAB was used as the chromogen.

Counter stained with hematoxylin.

Positive/negative staining were presented.

Only the secondary antibody was used as the negative control.