

Anti-PCSK9 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P105088

Product Information


Application	WB, IF (Cell)/ICC, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IF 1:50~1:200
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-694 of mouse PCSK9. (NP_705793.1).
Format	Affinity purified polyclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PCSK9 Rabbit pAb is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	FH3; PC9; Narc1; HCHOLA3; PCSK9.
Calculated MW	Calculated MW: 74 kDa
Uniprot ID	Q80W65
Gene ID	100102
Background	Enables apolipoprotein binding activity; lipoprotein particle binding activity; and low-density lipoprotein particle receptor binding activity. Involved in several processes, including cellular response to insulin stimulus; cellular response to starvation; and regulation of neuron apoptotic process. Acts upstream of or within several processes, including low-density lipoprotein receptor particle metabolic process; regulation of low-density lipoprotein particle receptor catabolic process; and triglyceride metabolic process. Located in COPII-coated ER to Golgi transport vesicle; endoplasmic reticulum; and extracellular space. Is expressed in several structures, including alimentary system; cerebellum; genitourinary system; liver; and telencephalon. Human ortholog(s) of this gene implicated in familial hypercholesterolemia and hypobetalipoproteinemia. Orthologous to human PCSK9 (proprotein convertase subtilisin/kexin type 9).



Immunofluorescence analysis of HeLa cells using PCSK9 Rabbit pAb(P105088) at a dilution of 1:200 (40× lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using PCSK9 Rabbit pAb(P105088) at a dilution of 1:200 (40× lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of lysates from Mouse liver, using PCSK9 Rabbit pAb (P105088) at 1:1,000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Kit.

Exposure time: 60s.