

Anti-Phospho-MAX (Ser11) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P108609

Product Information

| | |
|-------------|---|
| Application | WB, ELISA |
| Reactivity | Human |
| Dilution | WB 1:500~1:2,000 |
| Host | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Label | Unconjugated |
| Immunogen | A synthetic phosphorylated peptide around S11 of human MAX (NP_002373.3). |
| Format | Affinity purified polyclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3. |
| Storage | Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage. |
| Precautions | Anti-Phospho-MAX (Ser11) Rabbit pAb is for research use only and not for use in diagnostic or therapeutic procedures. |

Protein Information

| | |
|---------------|--|
| Synonyms | bHLHd4; Phospho-MAX-S11. |
| Calculated MW | Calculated MW: 18 kDa; Observed MW: 21 kDa |
| Uniprot ID | P61244 |
| Gene ID | 4149 |
| Background | The protein encoded by this gene is a member of the basic helix-loop-helix leucine zipper (bHLHZ) family of transcription factors. It is able to form homodimers and heterodimers with other family members, which include Mad, Mx1 and Myc. Myc is an oncoprotein implicated in cell proliferation, differentiation and apoptosis. The homodimers and heterodimers compete for a common DNA target site (the E box) and rearrangement among these dimer forms provides a complex system of transcriptional regulation. Mutations of this gene have been reported to be associated with hereditary pheochromocytoma. A pseudogene of this gene is located on the long arm of chromosome 7. Alternative splicing results in multiple transcript variants. |



Western blot analysis of lysates from Jurkat cells, using Phospho-MAX-S11 Rabbit pAb (P108609) at 1:1,000 dilution. Jurkat cells were treated by serum-starvation overnight. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% BSA.