

## Anti-IL-18 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P100682

### Product Information


Application	WB, IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IP 0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 37-193 of human IL18 (NP_001553.1).
Format	Affinity purified polyclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-IL-18 Rabbit pAb is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	IGIF; IL-18; IL-1g; IL1F4; IL18.
Calculated MW	Calculated MW: 22 kDa; Observed MW: 18 kDa
Uniprot ID	Q14116
Gene ID	3606
Background	The protein encoded by this gene is a proinflammatory cytokine of the IL-1 family that is constitutively found as a precursor within the cytoplasm of a variety of cells including macrophages and keratinocytes. The inactive IL-18 precursor is processed to its active form by caspase-1, and is capable of stimulating interferon gamma production, and of regulating both T helper (Th) 1 and Th2 responses. This cytokine has been implicated in the injury of different organs, and in potentially fatal conditions characterized by a cytokine storm. In humans, IL-18 gene is located on chromosome 11. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.



Western blot analysis of Recombinant Human IL-18 Protein using IL18 Rabbit pAb (P100682) at 1:500 dilution.  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.  
Lysates/proteins: 10 ng/1 ng per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Kit.  
Exposure time: 20s.



Immunoprecipitation analysis of 200 µg extracts of HeLa cells, using 3 µg IL18 antibody (P100682). Western blot was performed from the immunoprecipitate using IL18 antibody (P100682) at a dilution of 1:1000.