

Anti-Phospho-GAB1 (Tyr307) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P108595

Product Information

Application	WB, ELISA
Reactivity	Mouse
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic phosphorylated peptide around Y307 of human GAB1 (NP_002030.2).
Format	Affinity purified polyclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-GAB1 (Tyr307) Rabbit pAb is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	DFNB26; Phospho-GAB1-Y307.
Calculated MW	Calculated MW: 77 kDa; Observed MW: 115 kDa
Uniprot ID	Q13480
Gene ID	2549
Background	The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.



Western blot analysis of lysates from NIH/3T3 cells using Phospho-GAB1-Y307 Rabbit pAb (P108595) at 1:1,000 dilution. NIH/3T3 cells were treated by IGF-1 (50 ng/ml) at 37°C for 30 minutes after serum-starvation overnight.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Kit (SQ201).

Exposure time: 60s.