

## Anti-FOSL1 Rabbit mAb

Purified Recombinant Rabbit Monoclonal Antibody


Catalog # R015967

### Product Information

Application	WB, IF (Cell)/ICC, ELISA
Reactivity	Human
Dilution	WB 1:1,000~1:3,000; IF 1:100~1:200
Host	Rabbit
Clonality	Monoclonal
Clone No.	32A72M92
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human FRA1
Format	Affinity purified monoclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-FOSL1 Rabbit mAb [32A72M92] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

Synonyms	FRA1; FOSL1; Fos-related antigen 1; FRA-1.
Calculated MW	Calculated MW: 29 kDa; Observed MW: 40 kDa
Uniprot ID	P15407
Gene ID	8061
Background	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]
Cellular Location	Nucleus



Western Blot - Anti-FOSL1 Rabbit mAb [32A72M92]

All lanes: R015967 at 1:3,000 dilution

Lane 1: HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates

Lane 2: HepG2 (Human hepatocarcinoma epithelial cell) whole cell lysates

Lane 3: HCT116 (Human colorectal carcinoma epithelial cell) whole cell lysates

Lane 4: A431 (Human epidermoid teratoma cell line) whole cell lysates

Lane 5: K562 (Human chronic myeloid leukemia cell) whole cell lysates


Lysates/proteins at 10 µg per lane.

Secondary antibody: Goat Anti-Rabbit IgG (H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution

Predicted band size: 29 kDa

Observed band size: 40 kDa

Developed using the ECL technique (Cat. No. SQ201).



Immunofluorescence - Anti-FOSL1 Rabbit mAb [32A72M92]

Sample: HeLa cells

The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.

Primary antibodies: R015967 at 1:100 dilution and  $\alpha$ -tubulin Mouse Monoclonal Antibody (Cat. No. LF209) at 1:100 dilution

Secondary antibodies: Goat anti-Rabbit (488) at 1:1,000 dilution (shown in green) and Goat anti-Mouse (555) at 1:1,000 dilution (shown in red)

Nuclei were stained with DAPI (shown in blue).